

**Hanita lenses A constant conversion table**

		A (SRK II)	A (SRK/T)	pACD	sf	a0	a1	a2
<b>Aspheric AF</b>								
<b>SeeLens AF</b>	IOL MASTER	119.53	119.1	5.58	1.8	1.375	0.4	0.1
	US BIOMETRY	119.1	118.76	5.37	1.58	1.15	0.4	0.1
<b>BunnyLens AF</b>	IOL MASTER	118.78	118.5	5.2	1.42	0.978	0.4	0.1
	US BIOMETRY	118.35	118.16	4.98	1.2	0.753	0.4	0.1

<b>Hydrophobic HP</b>								
<b>SeeLens HP</b>	IOL MASTER	120.03	119.5	5.84	2.06	1.64	0.4	0.1
	US BIOMETRY	119.6	119.16	5.62	1.84	1.415	0.4	0.1

<b>Multifocal MF</b>								
<b>SeeLens MF</b>	IOL MASTER	118.9	118.6	5.26	1.48	1.044	0.4	0.1
	US BIOMETRY	118.48	118.26	5.05	1.27	0.819	0.4	0.1

<b>Spherical</b>								
<b>B-Lens</b>	IOL MASTER	118.83	118.54	5.23	1.44	1.004	0.4	0.1
	US BIOMETRY	118.4	118.2	5.01	1.23	0.779	0.4	0.1
<b>SeeLens</b>	IOL MASTER	118.9	118.6	5.26	1.48	1.044	0.4	0.1
	US BIOMETRY	118.48	118.26	5.05	1.27	0.819	0.4	0.1
<b>BunnyLens</b>	IOL MASTER	118.83	118.54	5.23	1.44	1.004	0.4	0.1
	US BIOMETRY	118.4	118.2	5.01	1.23	0.779	0.4	0.1